

Configuration Trieur SONY SH800 :

Lasers	PMT	Miroir	Filtre BP	Option de fluorochromes
Bleu 488 nm	FL1	NONE	525/50	EGFP, CFSE, FITC, Alexa Fluor 488, EYFP, mCitrine
	FL2	561LP	585/30	PE
	FL3	600LP	617/30	PE-Texas Red, PI
	FL4	639LP	665/30	7-AAD, PE-Cy5, PerCP
	FL5	685LP	720/60	PerCP-Cy5.5, PE-Cy5.5, PerCP-eFluor 710
	FL6	752LP	785/60	PE-Cy7
Jaune/Vert 561nm	FL2	561LP	585/30	mOrange, DsRed-Monomer, tdTomato,
	FL3	600LP	617/30	mCherry, PE-CF594
	FL4	639LP	665/30	PE-Texas Red, PI, mPlum